improving people's lives



September 5, 2013

Re: Sunrise Medical Power Wheelchair Articulating/Elevating Legrest (Option code 172H125 and 195H125)

To Whom It May Concern:

The purpose of this letter is to provide details regarding the center mount articulating/elevating platform legrest that Sunrise Medical makes available as an option on its power wheelchairs. Our center mount legrest allows the articulation function and the elevation function to operate independently. The independent functions address medical needs which are common among people with permanent disabilities.

In order for the function to operate independently there must be 2 actuators; one to operate the articulation feature and one to operate the elevation. This must be taken into consideration when configuring the power wheelchair and determining whether an expandable controller and drive thru are required.

Please let me know if you have any questions regarding this information.

Regards,

Rita Hostak

Vice President Government Relations

ta A Hostak

Sunrise Medical