improving people's lives



September 5, 2013

Re: Sunrise Medical Power Wheelchair Recline with Power Shear Reduction (Option codes 172PA6, 172PA7, 172PA16, 172PA17, 105PT3, 105PT4)

To Whom It May Concern:

The purpose of this letter is to provide details regarding the Power Reclining Back with Power Shear Reduction that Sunrise Medical makes available as an option on its power wheelchairs. Our power reclining back with power shear reduction allows the shear reduction feature to operate independently from the reclining back.

In order for the functions to operate independently there must be 2 actuators; one to operate the reclining feature and one to operate the shear reduction feature. This must be taken into consideration when configuring the power wheelchair and determining whether an expandable controller and drive thru are required.

Please let me know if you have any questions regarding this information.

Regards,

Rita Hostak

Vice President Government Relations

Ita S. Hostak

Sunrise Medical