WHEELCHAIR SMARTER.

AWARD-WINNING SMART WHEELCHAIR TECHNOLOGY REIMAGINING MOBILITY AT EVERY TURN.

LUCI

LUCI is the first and only smart technology accessory for power wheelchairs designed to give riders a safer and more inclusive experience.

LUCI prevents power wheelchair tips, collisions, and falls before they happen, thanks to a platform which helps the wheelchair "see" its environment and respond to steering input from the user. The result is safe, confident independence – not autonomy.

AWARDED BY









LEARN MORE ABOUT LUCI + QUICKIE


SMART WHEELCHAIR TECHNOLOGY FOR THE MOST FEARLESS PEOPLE ON EARTH.



FUELED BY COLLABORATION

We are founded and led by committed members of the community of people living with disabilities and their caregivers. Our team is motivated by a singular mission to make this beautiful, dangerous world more accessible for all people using power wheelchairs. Along the way, we have created a culture of intentional partnership because technology is another area where inclusion matters.

LUCI and the MyLUCI app are compatible with:



LUCI'S 360° VIEW

LUCI adds hardware and software to riders' new and existing chairs, offering unparalleled stability, security, and connectivity.

Our patented system combines data from sensor fusion, infrared, ultrasonics, and radar into a single view of the world, enabling never-before-seen possibilities for power wheelchair riders.

LUCI is now available on top-selling models of QUICKIE and ZIPPIE power wheelchairs. For veterans within the Veterans Health Administration (VA), LUCI can be installed on new QUICKIE power wheelchairs as they are being assembled at the Sunrise Medical production facility.

Q700 M | Q500 M | Q300 M/ZIPPIE

