



JAY[®] Zip Pediatric Backrest

E2613

January 2025

Mark For:

Submitting for:

Quote Order

Date: _____
 Dealer Acct #: _____
 Dealer: _____
 Dealer Contact: _____
 Dealer Address: _____
 Dealer City: _____ ST: _____ ZIP: _____
 Dealer Phone: () Fax: ()
 Confirmation Email: _____
 Confirm Via: Fax Email

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____
 Attention: _____
 Address: _____
 Address: _____
 Ship To City: _____ ST: _____ ZIP: _____
 Ship To Phone: () Fax: ()

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

JAY ZIP BACK - E2613

<input type="checkbox"/>	JZIPBXXY	E2613	JAY Zip Backrest All Sizes	\$589.00
<i>NOTE: Select appropriate assembly number below. 165 lbs Weight Capacity</i>				

EACH BACK SHIPS WITH THE FOLLOWING:

- Ultra lightweight mounting hardware
- 2" of anterior depth adjustment, 2.75" with reversing mounting pins
- 10mm Combo Wrench
- 4mm Hex Key
- Neoprene Backrest Hardware Covers (Allows for quick release)
- Water resistant inner cover & Ionic+™ outer cover
- Adjustable lumbar support



BACKREST WIDTH	8"	10"	12"	14"
FITS BACK CANE WIDTHS	8"-9"	10"-11"	12"-13"	14"-15"
MAX. USER TRUNK WIDTH	7"	9"	11"	13"
ACTUAL SHELL WIDTH	5.6"	7.6"	9.6"	11.6"
DEPTH OF CONTOUR	1"	1.5"	2"	2.5"

		Backrest Width			
		8"	10"	12"	14"
Backrest Height	6" (S)	<input type="checkbox"/> JZIPB08S	<input type="checkbox"/> JZIPB10S	<input type="checkbox"/> JZIPB12S	<input type="checkbox"/> JZIPB14S
	8" (M)	<input type="checkbox"/> JZIPB08M	<input type="checkbox"/> JZIPB10M	<input type="checkbox"/> JZIPB12M	<input type="checkbox"/> JZIPB14M
	10" (T)	<input type="checkbox"/> JZIPB08T	<input type="checkbox"/> JZIPB10T	<input type="checkbox"/> JZIPB12T	<input type="checkbox"/> JZIPB14T
	12" (X)	<input type="checkbox"/> JZIPB08X	<input type="checkbox"/> JZIPB10X	<input type="checkbox"/> JZIPB12X	<input type="checkbox"/> JZIPB14X

The JAY ZIP Backrest has been dynamically tested and meets all criteria for use in motor vehicles according to Section 5.1 of ISO 16840-4 and Section 20, ANSI/RESNA WC/Volume 4.
 Please visit www.sunparts.us for spare parts information.

JAY ZIP BACK ACCESSORIES

JAY ZIP Back accessories are optional. All accessories will be factory mounted as much as possible to ensure proper shipping.

LATERAL SUPPORTS

Use this section for lateral supports.

FIXED LATERAL THORACIC SUPPORT

Includes fixed mounting hardware.


<input type="checkbox"/>	JZFX34F	E0956	LAT ZIP FIXED 3X4 FLAT	\$292.00	<input type="checkbox"/>	JZFX34C	E0956	LAT ZIP FIXED 3X4 CURVE	\$303.00
<input type="checkbox"/>	JZFX45F	E0956	LAT ZIP FIXED 4X5 FLAT	\$292.00	<input type="checkbox"/>	JZFX45C	E0956	LAT ZIP FIXED 4X5 CURVE	\$303.00

SWING-AWAY LATERAL THORACIC SUPPORT

Includes swing-away mounting hardware and all laterals are sold in pairs.

<input type="checkbox"/>	JZSA34F		LAT 3X4 FLAT W/ S/A MOUNT	\$904.00	<input type="checkbox"/>	JZSA34CO		LAT 3X4 CRV 1" OFF W/ S/A MNT	\$915.00
	JZ34F	E0956	LAT ZIP 3X4 FLAT	\$292.00		JZ34CO	E0956	LAT 3X4 CURVE W/ 1" OFFSET	\$303.00
	JZSAMNT	E1028	LATERAL SWING AWAY MOUNT	\$612.00		JZSAMNT	E1028	LATERAL SWING AWAY MOUNT	\$612.00
<input type="checkbox"/>	JZSA45F		LAT 4X5 FLAT W/ S/A MOUNT	\$904.00	<input type="checkbox"/>	JZSA45CO		LAT 4X5 CRV 1" OFF W/ S/A MNT	\$915.00
	JZ45F	E0956	LAT ZIP 4X5 FLAT	\$292.00		JZ45CO	E0956	LAT 4X5 CURVE W/ 1IN OFFSET	\$303.00
	JZSAMNT	E1028	LATERAL SWING AWAY MOUNT	\$612.00		JZSAMNT	E1028	LATERAL SWING AWAY MOUNT	\$612.00
<input type="checkbox"/>	JZSA34C		LAT 3X4 CURVE W/ S/A MOUNT	\$915.00	<input type="checkbox"/>	JZSA34FOA		LAT 3X4 FLAT ADJ 1" OFF W/ S/A	\$910.00
	JZ34C	E0956	LAT ZIP 3X4 CURVE	\$303.00		JZ34FOA	E0956	LAT 3X4 FLAT W/ ADJ 1" OFFSET	\$298.00
	JZSAMNT	E1028	LATERAL SWING AWAY MOUNT	\$612.00		JZSAMNT	E1028	LATERAL SWING AWAY MOUNT	\$612.00
<input type="checkbox"/>	JZSA45C		LAT 4X5 CURVE W/ S/A MOUNT	\$915.00	<input type="checkbox"/>	JZSA45FOA		LAT 4X5 FLAT ADJ 1" OFF W/ S/A	\$910.00
	JZ45C	E0956	LAT ZIP 4X5 CURVE	\$303.00		JZ45FOA	E0956	LAT 4X5 FLAT W/ ADJ 1" OFFSET	\$298.00
	JZSAMNT	E1028	LATERAL SWING AWAY MOUNT	\$612.00		JZSAMNT	E1028	LATERAL SWING AWAY MOUNT	\$612.00
<input type="checkbox"/>	JZSA34FO		LAT 3X4 FLT 1" OFF W/ S/A MNT	\$904.00	<input type="checkbox"/>	JZSA34COA		LAT 3X4 CRV ADJ 1" OFF W/ S/A	\$921.00
	JZ34FO	E0956	LAT ZIP 3X4 FLAT W/ 1" OFFSET	\$292.00		JZ34COA	E0956	LAT 3X4 CRV W/ ADJ 1" OFFSET	\$309.00
	JZSAMNT	E1028	LATERAL SWING AWAY MOUNT	\$612.00		JZSAMNT	E1028	LATERAL SWING AWAY MOUNT	\$612.00
<input type="checkbox"/>	JZSA45FO		LAT 4X5 FLT 1" OFF W/ S/A MNT	\$904.00	<input type="checkbox"/>	JZSA45COA		LAT 4X5 CRV ADJ 1" OFF W/ S/A	\$921.00
	JZ45FO	E0956	LAT ZIP 4X5 FLAT W/ 1" OFFSET	\$292.00		JZ45COA	E0956	LAT 4X5 CRV W/ ADJ 1" OFFSET	\$309.00
	JZSAMNT	E1028	LATERAL SWING AWAY MOUNT	\$612.00		JZSAMNT	E1028	LATERAL SWING AWAY MOUNT	\$612.00


LATERAL MEASURING GUIDE (Fixed and Swing Away)



		Back Width			
DIM	Description	8" Wide	10" Wide	12" Wide	14" Wide
A	Lateral Adjustment Range (Min Spacing - Max Spacing)	3.8" - 4.8"	5.1" - 6.1"	6.8" - 7.8"	8.4" - 9.4"
B	Lateral Height Range (Measured from bottom edge of shell to top of lateral)	6" Tall Back			
		8" Tall Back			
		10" Tall Back			
		12" Tall Back			

HEADRESTS

<input type="checkbox"/>	JZHED6		Headrest w/ Pad - 6" Plush	\$343.00
	JPLSH06	E0955	Pad - 6" Plush	\$153.00
	JZIPBHH	E1028	Onyx post w/ Axy's Hardware	\$190.00
<input type="checkbox"/>	JZHED8		Headrest w/ Pad - 8" Plush	\$343.00
	JPLSH08	E0955	Pad - 8" Plush	\$153.00
	JZIPBHH	E1028	Onyx post w/ Axy's Hardware	\$190.00
<input type="checkbox"/>	JZHEDCINF		Headrest w/ Pad - 6" Contoured	\$326.00
	JCNTR06	E0955	Pad - 6" Contoured Cradle	\$136.00
	JZIPBHH	E1028	Onyx post w/ Axy's Hardware	\$190.00
<input type="checkbox"/>	JZHEDCS		Headrest w/ Pad - 8" Contoured	\$337.00
	JCNTR08	E0955	Pad - 8" Contoured Cradle	\$147.00
	JZIPBHH	E1028	Onyx post w/ Axy's Hardware	\$190.00



ONYX Post w/ AXYS Hdw & Adptr



Plush Pad



Contoured Pad

Vanity Flap







Vanity flaps are NOT INCLUDED as a standard option with the JAY ZIP Back.

<input type="checkbox"/>	VF08	K0108	8" to 9" Wide	\$52.00	<input type="checkbox"/>	VF12	K0108	12" to 13" Wide	\$52.00
<input type="checkbox"/>	VF10	K0108	10" to 11" Wide	\$52.00	<input type="checkbox"/>	VF14	K0108	14" to 15" Wide	\$52.00

Spine Align Positioning Componets

<input type="checkbox"/>	J3SPNCR	K0108	Curved Shim	\$29.00
<input type="checkbox"/>	J3SPNBS	K0108	Small Bean Bag	\$32.00
<input type="checkbox"/>	J3SPNFL	K0108	Flat Shim	\$29.00

<input type="checkbox"/>	J3SPNBL	K0108	Large Bean Bag	\$35.00
<input type="checkbox"/>	J3SPNLS	K0108	Lateral Shim	\$55.00
<input type="checkbox"/>	J3SPNLP	K0108	Posterior/Lateral Pelvic	\$55.00




						
	Curved Shim	Flat	Lateral Shim	Sm Bean Bag	Lg Bean Bag	Post/Lat Pelvic
W	10"	10"	2.4"	6"	10"	5.3"
H	4"	4"	4.3"	4.5"	3.2"	3.6"
D	0.5"	0.3"	0.5"	0.75"	0.75"	0.75"

CHEST STRAPS & HARNESSSES

Chest Straps & Harnesses include guide brackets and hardware.

<input type="checkbox"/>	JZ9407	E0960	Center Opening w/ Buckle (EI)	\$177.00
<input type="checkbox"/>	JZ9400	E0960	Center Opening w/ Buckle (XS)	\$177.00
<input type="checkbox"/>	JZ9401	E0960	Center Opening w/ Buckle (S)	\$177.00
<input type="checkbox"/>	JZ9402	E0960	Center Opening w/ Buckle (M)	\$177.00
<input type="checkbox"/>	JZ8107	E0960	Classic - Dynamic Pull Down (EI)	\$177.00
<input type="checkbox"/>	JZ8100	E0960	Classic - Dynamic Pull Down (XS)	\$177.00
<input type="checkbox"/>	JZ8101	E0960	Classic - Dynamic Pull Down (S)	\$177.00
<input type="checkbox"/>	JZ8102	E0960	Classic - Dynamic Pull Down (M)	\$177.00
<input type="checkbox"/>	JZ8300	E0960	Classic - Structured Pull Down (XS)	\$177.00
<input type="checkbox"/>	JZ8301	E0960	Classic - Structured Pull Down (S)	\$177.00
<input type="checkbox"/>	JZ8302	E0960	Classic - Structured Pull Down (M)	\$177.00

<input type="checkbox"/>	JZ9107	E0960	Contoured - Dynamic Pull Down (EI)	\$177.00
<input type="checkbox"/>	JZ9100	E0960	Contoured - Dynamic Pull Down (XS)	\$177.00
<input type="checkbox"/>	JZ9101	E0960	Contoured - Dynamic Pull Down (S)	\$177.00
<input type="checkbox"/>	JZ9102	E0960	Contoured - Dynamic Pull Down (M)	\$177.00
<input type="checkbox"/>	JZ9300	E0960	Contoured - Strctrd Pull Down (XS)	\$177.00
<input type="checkbox"/>	JZ9301	E0960	Contoured - Strctrd Pull Down (S)	\$177.00
<input type="checkbox"/>	JZ9302	E0960	Contoured - Strctrd Pull Down (M)	\$177.00
<input type="checkbox"/>	JZ4050	E0960	Padded Anterior Trunk Straps (XS)	\$165.00
<input type="checkbox"/>	JZ4000	E0960	Padded Anterior Trunk Straps (S)	\$165.00
<input type="checkbox"/>	JZ4100	E0960	Padded Anterior Trunk Straps (M)	\$165.00

																																																																																																																
Classic Style	Contour Style	Center-Opening with Buckle																																																																																																														
<table border="1"> <thead> <tr> <th></th> <th>A</th> <th>B</th> <th>C</th> <th>D</th> <th>E</th> <th>F</th> <th>G</th> <th>H</th> </tr> </thead> <tbody> <tr> <td>EI</td> <td>9.8"</td> <td>3"</td> <td>2"</td> <td>5.3"</td> <td>3"</td> <td>3"</td> <td>1"</td> <td>6.5"</td> </tr> <tr> <td>XS</td> <td>12.5"</td> <td>3.5"</td> <td>2"</td> <td>7.5"</td> <td>3.5"</td> <td>3"</td> <td>2"</td> <td>7.5"</td> </tr> <tr> <td>S</td> <td>15"</td> <td>4.5"</td> <td>2"</td> <td>9"</td> <td>4"</td> <td>4"</td> <td>2"</td> <td>8.5"</td> </tr> <tr> <td>M</td> <td>19"</td> <td>5.5"</td> <td>2"</td> <td>11"</td> <td>5"</td> <td>5"</td> <td>3"</td> <td>9.5"</td> </tr> </tbody> </table>		A	B	C	D	E	F	G	H	EI	9.8"	3"	2"	5.3"	3"	3"	1"	6.5"	XS	12.5"	3.5"	2"	7.5"	3.5"	3"	2"	7.5"	S	15"	4.5"	2"	9"	4"	4"	2"	8.5"	M	19"	5.5"	2"	11"	5"	5"	3"	9.5"	<table border="1"> <thead> <tr> <th></th> <th>A</th> <th>B</th> <th>C</th> <th>D</th> <th>E</th> <th>F</th> <th>G</th> <th>H</th> </tr> </thead> <tbody> <tr> <td>EI</td> <td>9.8"</td> <td>3"</td> <td>1.8"</td> <td>5.3"</td> <td>1"</td> <td>3.5"</td> <td>1.5"</td> <td>6.5"</td> </tr> <tr> <td>XS</td> <td>12.5"</td> <td>3.5"</td> <td>2"</td> <td>7.5"</td> <td>1"</td> <td>3.5"</td> <td>1.5"</td> <td>7.5"</td> </tr> <tr> <td>S</td> <td>15"</td> <td>4.5"</td> <td>2"</td> <td>9"</td> <td>1"</td> <td>4"</td> <td>2"</td> <td>8.5"</td> </tr> <tr> <td>M</td> <td>18.5"</td> <td>5.5"</td> <td>2"</td> <td>11"</td> <td>1"</td> <td>5"</td> <td>2.5"</td> <td>9.5"</td> </tr> </tbody> </table>		A	B	C	D	E	F	G	H	EI	9.8"	3"	1.8"	5.3"	1"	3.5"	1.5"	6.5"	XS	12.5"	3.5"	2"	7.5"	1"	3.5"	1.5"	7.5"	S	15"	4.5"	2"	9"	1"	4"	2"	8.5"	M	18.5"	5.5"	2"	11"	1"	5"	2.5"	9.5"	<table border="1"> <thead> <tr> <th></th> <th>A</th> <th>B</th> <th>C</th> </tr> </thead> <tbody> <tr> <td>EI</td> <td>9.8"</td> <td>7"</td> <td>5"</td> </tr> <tr> <td>XS</td> <td>11"</td> <td>8"</td> <td>5.8"</td> </tr> <tr> <td>S</td> <td>12.3"</td> <td>9"</td> <td>6.5"</td> </tr> <tr> <td>M</td> <td>14.3"</td> <td>11"</td> <td>9.8"</td> </tr> </tbody> </table>		A	B	C	EI	9.8"	7"	5"	XS	11"	8"	5.8"	S	12.3"	9"	6.5"	M	14.3"	11"	9.8"
	A	B	C	D	E	F	G	H																																																																																																								
EI	9.8"	3"	2"	5.3"	3"	3"	1"	6.5"																																																																																																								
XS	12.5"	3.5"	2"	7.5"	3.5"	3"	2"	7.5"																																																																																																								
S	15"	4.5"	2"	9"	4"	4"	2"	8.5"																																																																																																								
M	19"	5.5"	2"	11"	5"	5"	3"	9.5"																																																																																																								
	A	B	C	D	E	F	G	H																																																																																																								
EI	9.8"	3"	1.8"	5.3"	1"	3.5"	1.5"	6.5"																																																																																																								
XS	12.5"	3.5"	2"	7.5"	1"	3.5"	1.5"	7.5"																																																																																																								
S	15"	4.5"	2"	9"	1"	4"	2"	8.5"																																																																																																								
M	18.5"	5.5"	2"	11"	1"	5"	2.5"	9.5"																																																																																																								
	A	B	C																																																																																																													
EI	9.8"	7"	5"																																																																																																													
XS	11"	8"	5.8"																																																																																																													
S	12.3"	9"	6.5"																																																																																																													
M	14.3"	11"	9.8"																																																																																																													

JAY YOUR WAY MODIFICATIONS

Note: To order options offered via JAY Your Way add an "M" to the front of the part number and follow the prompts.

For additional modifications visit www.sunrisemedical.com to view even more!

COVER COLOR OPTIONS - Rear Center Flap

<input type="checkbox"/>	COCPI	Pink Semi Stretch Cover	NC
<input type="checkbox"/>	COCB	Blue Semi Stretch Cover	NC
<input type="checkbox"/>	COCPU	Purple Cover	NC
<input type="checkbox"/>	COCR	Red Cover	NC
<input type="checkbox"/>	COCB	Toxic Green Cover	NC

EMBROIDERY

<input type="checkbox"/>	EMBROIDERY	K0108	Embroidery (32 Character Max)	\$26.00
--------------------------	------------	-------	-------------------------------	---------

Location: Front
 Font: Block Script

Embroidery Text:

COLOR BACK SHELLS

Paint your Zip back shell to match any of the Zippie or Quickie Frame Colors!

<input type="checkbox"/>	XXCBS	Zippie/Quickie Back Shell Paint	NC
--------------------------	-------	---------------------------------	----

Write Desired Back Shell Color Here:

SPINAL FLUID PAD

<input type="checkbox"/>	SFPS	Small Spinal Fluid Pad	\$58.00
<input type="checkbox"/>	SFPM	Medium Spinal Fluid Pad	\$58.00

*Don't see the option you want? NO PROBLEM!
 Write your request below and we can evaluate it and provide a quote.*

Notes:



Sunrise Medical (US) LLC
 12002 Volunteer Blvd. · Mount Juliet, TN 37122 · USA