Quickie SR45™
May 2025
EIZ9

Mark For:

Date:
Dealer Acct #:
Dealer:
Dealer Contact:
Dealer Address:
Dealer City: ST: ZIP:
Dealer Phone: ( ) Fax: ( )
Confirmation Email:
Confirm Via: ☐ Fax ☐ Email

Submitting for: ☐ **Quote** ☐ **Order**

PO#:

ADDITIONAL SHIPPING INFORMATION

Ship To:
Attention:
Address:
Address:
Ship To City: ST: ZIP:
Ship To Phone: ( ) Fax: ( )

*The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.*

1. FRAME

Base Model

✔	Option #	HCPCS	Definition	Price
☐	194WM77	E1161	Quickie SR45 *For 14" sizes use E1233*	$3,198.00

Power Tilt

✔	Option #	HCPCS	Definition	Price
☐	194PA53	K0108 E1002	Power Tilt *(W/ Attendant &Tilt Cntrl)* Tilt Control mounts to armrest ☐ Left ☐ Right (Rear Facing)	$2,964.00

Frame Options

✔	Option #	HCPCS	Definition	Price
☐	194FO3		Standard Base and Seat Frame *Maximum Weight: 265 lbs.*	NC

Transit

✔	Option #	HCPCS	Definition	Price
☐	194TO38	K0108	Transit Option	$374.00

Tilt Release

✔	Option #	HCPCS	Definition	Price
☐	194TT12		Backrest Dual Trigger Cable Release *(N/A w/ Power Tilt)*	NC
☐	194TT13		Foot Tilt Release - *Mounted on Right Side of frame (N/A w/ Power Tilt)*	$58.00

Tilt Angle

✔	Option #	HCPCS	Definition	Price
☐	194TA2		0° to 45°	NC



Quickie SR45 shown with the JAY© J3™ backrest and JAY© 3™ cushion.



Quickie SR45 features IRIS™ (Intelligent Rotation In-Space) technology.

© 5/2025 MK-100025 REV Q

✔	Option #	HCPCS	Definition	Price	✔	Option #	HCPCS	Definition	Price

2. SIZE / SEATING

Frame Width

☐	194W5		14"	NC
☐	194W6		15"	NC
☐	194W7		16"	NC
☐	194W8		17"	NC
☐	194W9		18"	NC
☐	194W10		19"	NC
☐	194W11	E2201	20"	$583.00

Seat Type

☐	194ST35	E2231	Standard Seat Pan	$234.00
☐	194ST7		Omit Seat Pan	NC

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list.
Cushion name or part number must be filled out with desired width and depth

		JAY Bundled Cushion		
☐	194JC20	Cushion Name or Part #		Price List
		Width:	Depth:	

Frame Depth Set Up (Fill out either A or B)

A. Setup By QUICKIE

A. Desired Seat Depth

See **Figure 1** for pricing, **Figure 2**, Dimension A for measurement

☐	194SD7		14"	NC
☐	194SD8		15"	NC
☐	194SD9		16"	NC
☐	194SD10		17"	Figure 1.1
☐	194SD11		18"	Figure 1.1
☐	194SD12		19"	Figure 1.1
☐	194SD13	E2203	20"	Figure 1.1
☐	194SD14	E2203	21"	Figure 1.1
☐	194SD33	E2204	22"	Figure 1.1

A. User's Seated Height

See **Figure 2**, Dimension B. Distance from seat surface to top of the head.

☐	194SJ01		Less than 37"	NC
☐	194SJ02		37" or Greater	NC

OR

B. Setup By DEPTH

B. Frame Depth Range

☐	194FL91		Short Seat Frame 14"-18"	NC
☐	194FL92		Medium Seat Frame 16"-20"	NC
☐	194FL93		Long Seat Frame 18"-22"	NC

B. Frame Depth

☐	194FD5		14" *recommend short seat frame*	NC
☐	194FD6		15" *recommend short seat frame*	NC
☐	194FD7		16" *recommend short seat frame*	NC
☐	194FD8		17" *recommend short seat frame*	NC
☐	194FD9		18" *recommend medium seat frame*	NC
☐	194FD10		19" *recommend medium seat frame*	NC
☐	194FD11	E2203	20" *recommend long seat frame*	$641.00
☐	194FD16	E2203	21" *recommend long seat frame*	$641.00
☐	194FD17	E2204	22" *recommend long seat frame*	$1,051.00

Figure 1

Desired Seat Depth	3	2	1	Flush	-1	-2
14"	NC	NC	NC	NC	N/A	N/A
15"	NC	NC	NC	NC	N/A	N/A
16"	NC	NC	NC	NC	NC	N/A
17"	$641	NC	NC	NC	NC	NC
18"	$641	$641	NC	NC	NC	NC
19"	$1,051	$641	$641	NC	NC	NC
20"	N/A	$1,051	$641	$641	NC	NC
21"	N/A	N/A	$1,051	$641	$641	NC
22"	N/A	N/A	N/A	$1,051	$641	$641

A. Top of Cushion to Seat Rail

See **FIGURE 3** Dimensions A & B. This measurement is taken from the top of the cushion to the upper surface of the seat rail. If the cushion has a wedged shape and/or has significant lateral pelvic contour, then determine the height of the cushion at the point where the ischial tuberosities would contact (not immerse) the cushion surface.

☐	194CJ04		2 1/2" or Greater	NC
☐	194CJ05		Less than 2 1/2"	NC
☐	194CJ03		Below Seat Rail	NC

Figure 2





Figure 3



© 5/2025 MK-100025 REV Q

✔	Option #	HCPCS	Definition	Price	✔	Option #	HCPCS	Definition	Price

3. CASTERS / SEAT HEIGHT

Rear Wheel Size	CASTER SIZE					
	☐ 4	☐ 5	☐ 6	☐ 8x1	☐ 8x1.5	☐ 8x2
☐ 12"		☐ 16.0	☐ 16.0			
		☐ 16.5	☐ 16.5			
		☐ 17.0	☐ 17.0			
		☐ 17.5	☐ 17.5	☐ 17.5	☐ 17.5	☐ 17.5
☐ 12" Low Poly		☐ 15.0				
		☐ 15.5	☐ 16.0			
		☐ 16.0	☐ 16.5			
		☐ 16.5	☐ 17.0			
☐ 16"	☐ 14.0	☐ 15.0				
	☐ 15.0	☐ 17.0	☐ 17.0			
		☐ 17.5	☐ 17.5	☐ 17.5	☐ 17.5	☐ 17.5
		☐ 18.0	☐ 18.0	☐ 18.0	☐ 18.0	
		☐ 18.5	☐ 18.5	☐ 18.5	☐ 18.5	☐ 18.5
					☐ 19.0	
☐ 20"		☐ 16.0	☐ 16.0			
		☐ 16.5	☐ 16.5			
		☐ 17.0	☐ 17.0		☐ 17.5	
		☐ 17.5	☐ 17.5	☐ 19.5	☐ 20.0	☐ 17.5
			☐ 19.5	☐ 20.0	☐ 20.5	☐ 19.5
☐ 22"		☐ 16.5	☐ 16.5			
		☐ 17.0	☐ 17.0		☐ 17.5	
		☐ 17.5	☐ 17.5	☐ 17.5	☐ 18.0	☐ 17.5
		☐ 18.0	☐ 18.0	☐ 18.0	☐ 18.5	
☐ 24"		☐ 17.5	☐ 17.5	☐ 17.5	☐ 18.0	
		☐ 18.0	☐ 18.0	☐ 18.0	☐ 18.5	☐ 17.5
		☐ 18.5	☐ 18.5	☐ 18.5	☐ 19.0	☐ 18.5
		☐ 19.0	☐ 19.0	☐ 19.0	☐ 19.5	

Casters

☐	194CW73		4" x 1.25" Semi Pneumatic	NC
☐	194CW5		5" Low Profile Polyurethane	NC
☐	194CW7	E2219	5" X 1.5" Semi Pneumatic	$100.00
☐	194CW34		6" Polyurethane	NC
☐	194CW35	E2219	6"x 1.5" Semi Pneumatic	$100.00
☐	194CW32		8" Polyurethane	NC
☐	194CW39	E2219	8"x 1.5" Semi Pneumatic	$100.00
☐	194CW36	E2214 E2215	8"x 2" Pneumatic	$169.00
☐	194CW100	E2214 K0108	8"x 2" Pneumatic w/Insert	$292.00

4. REAR WHEELS / SEAT HEIGHT

Rear Wheel

☐	194RW13		12" Mag Wheel	NC
☐	194RW14		16" Mag Wheel	NC
☐	194RW3		20" Mag Wheel	NC
☐	194RW8		22" Mag Wheel	NC
☐	194RW1		24" Mag Wheel	NC

Rear Wheel Tire

☐	194RT3		Full Polyurethane *N/A with 12" and 16"*	NC
☐	194RT4		Low Polyurethane *Only available with 12" Mag Wheel*	NC
☐	194RT1		Pneumatic	NC
☐	194RT2	E2213 E2211	Pneumatic Airless Insert	$100.00

Handrim

12" and 16" wheels not available with handrims.

☐	194HR1		Aluminum Anodized	NC
☐	194HR4		Plastic Coated	$129.00
☐	194HR7		Omit Handrim	NC

Axle

☐	194AX1		Quick Release	NC
☐	194AX44		Threaded	$30.00

5. WHEEL LOCKS

Wheel Lock

☐	194WL8		Push To Lock	NC
☐	194WL43		Rear Foot Operated Wheel Lock	$88.00
☐	194WL10		Omit Wheel Lock	NC

Wheel Lock Extension Handle

Not available with rear foot operated wheel lock

☐	194EH7	E0961	Extension Handle (Pair)	$88.00

© 5/2025 MK-100025 REV Q

✔	Option #	HCPCS	Definition	Price	✔	Option #	HCPCS	Definition	Price

6. HANGER AND FOOTREST

Hanger

Front mount hangers only available w/ Alum angle adj footplates

☐	194H139		70° Swing In/Swing Out -Ext. Mount	NC
☐	194H140		70° Swing In/Swing Out -Front Mount	NC
☐	194H141		80° Swing In/Swing Out -Ext.Mount *N/A with all 8" casters*	NC
☐	194H142		80° Swing In/Swing Out -Front Mount	NC
☐	194H159	K0053	Articulating Hanger *N/A w/ Platform Flip Up and Alum Angle Adj Front Mount footplates.*	$431.00
☐	194H160	E0990	Elevating Hanger - Ext Mount *N/A w/ Platform Flip Up and Alum Angle Adj Front Mount footplates.*	$367.00
☐	194H161	E0990	Elevating Hanger - Front Mount *Only available w/ front mount alum angle adj footplates and 14"-22" width*	$367.00
☐	194H144		Omit Hanger	NC

Articulating / Elevating Hanger Cover

☐	194DI47		Neoprene Cover	$397.00
☐	194DI48		Gel Pad Cover	$467.00

Extension Tube (NA w/ Front Mount Hangers)

See lower leg length matrix

☐	194E1		Short (N/A with ELR)	NC
☐	194E2		Medium	NC
☐	194E5		Long	NC

Lower Leg Length Range Matrix

Hanger	Ext Tube	Extension Mount / Footplate: Composite	Alum Angle Adj	Aluminum
70° Swing In/Out	Short	10.0" - 14.0"	10.2" - 14.1"	NA
	Medium	11.0" - 15.9"	11.1" - 16.0"	NA
	Long	12.9" - 17.8"	13.0" - 17.9"	NA
80° Swing In/Out	Short	9.9" - 11.0"	10.1" - 13.7"	NA
	Medium	11.0" - 15.5"	11.1" - 15.7"	NA
	Long	12.7" - 17.3"	12.9" - 17.4"	NA

All values are ± 0.25" and based on flat seat heights and do not account for the minimum 2" ground clearance or footplate thickness. Full range not available with all seat heights. ELR measurements taken with ELR fully

Front Mount Extension Tubes (NA w/ Extension Mount Hangers)

Optional, not required for Front Mount Hangers

☐	194E23		2"	$63.00
☐	194E24		4" *N/A w/14", 15" Seat To Floor Ht, N/A w/80° Hangers & 8" casters*	$63.00

Hanger Mount			
Ext	**None**	**2"**	**4"**
70°	3.1" - 8.8"	3.1" - 11.1"	3.1" - 13.0"
80°	3.3" - 9.2"	3.3" - 11.7"	3.3" - 13.7"
90°	2.7" - 8.7"	2.7" - 11.2"	2.7" - 13.2"

Hanger Mounts

Possible interference with 80° & 90° hangers

☐	194MZ02		Standard Hanger Mount	STD
☐	194MZ01		1" Offset Hanger Mount	$36.00

Footplates

Heel loop not included, see Footplate Options section below

☐	194F1		Composite	NC
☐	194F5		Aluminum *Only Avail w/ Elevating Legrest*	NC
☐	194F291	K0040	Aluminium Angle Adjustable	$221.00

Footplate Options

☐	194HL1	E0951	Heel Loops	$47.00
☐	1945A02	E0952	Toe Loops (pair) *Only Avail w/ Angle Adj. Footplate*	$23.00

Hanger	Extension	Seat Sling to Top of Footrest: Min	Max	Min Allowable FSH
ALR	Short	8 1/2	16	10 1/2
	Med	10	18	12
	Long	12	20	14
ELR Ext. Mount	Short	8	14 1/2	10
	Med	8	16 1/2	10
	Long	9 1/2	18 1/2	11 1/2
ELR Front Mount	No Ext (top hole)	3 1/2	7 1/2	5 1/2
	2" ext tube	7 1/4	10 1/4	9 1/4
	4" ext tube	9 1/4	12 1/4	11 1/4
	6" ext tube	11 1/4	14 1/4	13 1/4

7. BACKREST / SEATING

Backrest Type

☐	194BT3		Non Folding Backrest *N/A with Height Adjustable Back*	NC
☐	194BT99		Folding Backrest *N/A with Reclining Backrest*	$58.00

Backrest Options

☐	194OD19		Removable Extension *N/A with Height adj. back with adjustable stroller handle.*	$100.00

JAY J3 Back

☐	194JB01		JAY J3 Back	See OF

Backpost Style - Select only one option within A, B, or C

A. Stroller Handle Backpost (Fixed)

Folding Backrest adds 1" to the listed values.

☐	194BP48		Low (18")	NC
☐	194BP49		Medium (21")	NC
☐	194BP50		Tall (24")	NC

OR

B. Height Adjustable Back w/ Adj. Stroller Handles


☐	194BP51		Standard 15"-21"	NC
☐	194BP52		Tall 18"-24"	NC

OR


C. Reclining Backrest

Not available with 22" Deep Frame

☐	194BP53	E1225	Medium (21")	$1,032.00
☐	194BP54	E1225	Tall (24")	$1,032.00



Reclining/Fixed Handle Ht. Diagram



Height Adj. Handle Height Diagram

© 5/2025 MK-100025 REV Q

✔	Option #	HCPCS	Definition	Price	✔	Option #	HCPCS	Definition	Price

7. BACKREST / SEATING (CONT)

A. Backrest Position (set up by Quickie)

See **FIGURE 4** Dimension C. The location of the backrest surface in relation to the front of the backpost. The measurement is taken from the front surface of the backpost to the deepest area of the backrest in contact with the user's body.

☐	194BQ01		+3"	NC
☐	194BQ05		+2"	NC
☐	194BQ02		+1"	NC
☐	194BQ03		Flush	NC
☐	194BQ04		-1"	NC
☐	194BQ06		-2"	NC




Figure 4

Back Type	
1	3" forward (+3")
2	2" forward (+2")
3	1" forward (+1")
4	Flush
5	1" rearward (-1")
6	2" rearward (-2")

8. ARMREST / SIDEGUARDS

Armrest

Armrests w/ Transfer Bar N/A w/ 14" Seat Depth

☐	194AA100	E0973	Dual Post Height Adj No Transfer Bar	$298.00
☐	194AA101	E0973	Dual Post Height Adj w/Transfer Bar	$309.00
☐	194AA102	E0973	Dual Post Flip Back Height Adj No Transfer Bar	$333.00
☐	194AA103	E0973	Dual Post Flip Back Height Adj w/Transfer Bar	$345.00
☐	194AA104	E0973	Single Post Height Adjustable - Std	$333.00
☐	194AA105	E0973	Single Post Height Adjustable - Low	$333.00
☐	194AA14		Omit Armrest	NC

Armrest Matrix for Single Post


Armrest Size	Transfer Tube Length D	Receiver	H Range	Pad Length L	
				Desk	Full
STD	8"	STD	10 3/4" to 15 3/4"	10 1/4"	14"
		LOW	8 3/4" to 13 3/4"		
LOW	5 1/4"	STD	9" to 13"		
		LOW	7 1/4" to 11 1/4"		

Armpad

☐	194AM25		Desk	NC
☐	194AM26		Full	NC

Armrest Receiver (Single Post Height Adjustable Only)

☐	194AR17		Low	NC
☐	194AR18		Std	NC



9. KEY OPTIONS AND ACCESSORIES

Anti Tips

☐	194AT1	E0971	Anti Tips (Pair)	$175.00

Positioning Belt

☐	194A141	E0978	1.5" Auto Buckle Universal w/Clamp	$82.00
☐	194A118	E0978	2" Aircraft Buckle Pad w/Clamp	$139.00
☐	194A120		Positioning Belt Clamps Only	$58.00

JAY Positioning

☐	194K1	E1020	JAY Lower Extremity S/A support ☐ Left ☐ Right Position	$460.00
☐	194K2		Omit Hanger ☐ Left ☐ Right	-$10.00
☐	194K7		Pad Size ☐ 6" x 8" ☐ 8" x 10"	NC

Misc.

☐	194A1		Black Utility Backpack	$88.00
☐	194TS03		Tilt Stops	$36.00
☐	194A999		Quickie Water Bottle & Holder	$63.00

© 5/2025 MK-100025 REV Q

✔	Option #	HCPCS	Definition	Price	✔	Option #	HCPCS	Definition	Price

10. COLOR

NEW!	Standard Frame Colors			
☐	194C214		Black (Gloss)	NC
☐	194C215		Black (Matte)	NC
☐	194C211		Celestial Sapphire	NC
☐	194C212		Desert Sage	NC
☐	194C213		Desert Tan	NC
☐	194C216		Metallic Shadow	NC
☐	194C217		Midnight Blue	NC
☐	194C218		Mystic Purple	NC
☐	194C219		Petal Pink	NC
☐	194C220		Pink Power	NC
☐	194C221		Red Pepper	NC
☐	194C222		Royal Blue	NC
☐	194C223		Salute Green	NC
☐	194C224		Scarlet Noir	NC
☐	194C225		Silver Ore	NC
☐	194C226		Solar Orange	NC
☐	194C227		Starlight White	NC
☐	194C228		Storm Gray	NC
☐	194C229		Turbo Green	NC
☐	194C230		Velocity Yellow	NC

Personalized Frame Color				
☐	194C112		Desert Camouflage	$460.00
☐	194C43		Mossy Oak Camouflage	$460.00
☐	194C42		Stars & Stripes	$460.00
☐	194C108		Zebra *(Choose above base color) Std with White unless above base paint is chosen N/A w/ Black, Matte Black)*	$460.00

Additional Options Through Built for Me

To inquire about any unique personalized Quickie products, contact: Phone: (800) 333 - 4000, Fax: (800) 333 9011

JAY® Signature Series Backrest and Accessories

Note: Not all sizes and configurations shown below, for full offering visit www.sunrisemedical.com


When ordering mark the part number desired and it will be added to the order as a separate line item .

JAY BACKREST

JAY J3 Posterior (PA) BACK - 2.2" Contour Depth

	Shell Height		Shell Width									
Mid Thoracic (MT)	12.0" (S)	J3PAMTSxxFS	☐ 12	$677	☐ 14	$677	☐ 16	$677	☐ 18	$677	☐ 20	$677
		J3PAMTSxxTH	☐ 12	$677	☐ 14	$677	☐ 16	$677	☐ 18	$677	☐ 20	$677
Upper Thoracic (UT)	16.5" (S)	J3PAUTSxxFS	☐ 12	$677	☐ 14	$677	☐ 16	$677	☐ 18	$677	☐ 20	$677
		J3PAUTSxxTH	☐ 12	$677	☐ 14	$677	☐ 16	$677	☐ 18	$677	☐ 20	$677
	18.0" (M)	J3PAUTMxxFS			☐ 14	$677	☐ 16	$677	☐ 18	$677	☐ 20	$677
		J3PAUTMxxTH			☐ 14	$677	☐ 16	$677	☐ 18	$677	☐ 20	$677
	19.5" (T)	J3PAUTTxxFS			☐ 14	$677	☐ 16	$677	☐ 18	$677	☐ 20	$677
		J3PAUTTxxTH			☐ 14	$677	☐ 16	$677	☐ 18	$677	☐ 20	$677
Shoulder Height (SH)	21.0" (S)	J3PASHSxxFS			☐ 14	$677	☐ 16	$677	☐ 18	$677	☐ 20	$677
		J3PASHSxxTH			☐ 14	$677	☐ 16	$677	☐ 18	$677	☐ 20	$677
	22.5" (M)	J3PASHMxxFS					☐ 16	$677	☐ 18	$677	☐ 20	$677
		J3PASHMxxTH					☐ 16	$677	☐ 18	$677	☐ 20	$677
	24.0" (T)	J3PASHTxxFS							☐ 18	$677	☐ 20	$677
		J3PASHTxxTH							☐ 18	$677	☐ 20	$677


FS Mounting Hardware



JAY J3 Posterior Lateral (PL) BACK - 3" Contour Depth

Mid Thoracic (MT)	15.0" (T)	J3PLMTTxxFS	☐ 12	$781	☐ 14	$781	☐ 16	$781	☐ 18	$781	☐ 20	$781
		J3PLMTTxxTH	☐ 12	$781	☐ 14	$781	☐ 16	$781	☐ 18	$781	☐ 20	$781
Upper Thoracic (UT)	16.5" (S)	J3PLUTSxxFS	☐ 12	$781	☐ 14	$781	☐ 16	$781	☐ 18	$781	☐ 20	$781
		J3PLUTSxxTH	☐ 12	$781	☐ 14	$781	☐ 16	$781	☐ 18	$781	☐ 20	$781
	18.0" (M)	J3PLUTMxxFS			☐ 14	$781	☐ 16	$781	☐ 18	$781	☐ 20	$781
		J3PLUTMxxTH			☐ 14	$781	☐ 16	$781	☐ 18	$781	☐ 20	$781
	19.5" (T)	J3PLUTTxxFS			☐ 14	$781	☐ 16	$781	☐ 18	$781	☐ 20	$781
		J3PLUTTxxTH			☐ 14	$781	☐ 16	$781	☐ 18	$781	☐ 20	$781
Shoulder Height (SH)	21.0" (S)	J3PLSHSxxFS			☐ 14	$781	☐ 16	$781	☐ 18	$781	☐ 20	$781
		J3PLSHSxxTH			☐ 14	$781	☐ 16	$781	☐ 18	$781	☐ 20	$781
	22.5" (M)	J3PLSHMxxFS					☐ 16	$781	☐ 18	$781	☐ 20	$781
		J3PLSHMxxTH					☐ 16	$781	☐ 18	$781	☐ 20	$781
	24.0" (T)	J3PLSHTxxFS							☐ 18	$781	☐ 20	$781
		J3PLSHTxxTH							☐ 18	$781	☐ 20	$781

TH Mounting Hardware



JAY J3 Posterior Deep (PD) BACK - 6" Contour Depth

Mid Thoracic (MT)	15.0" (T)	J3PDMTTxxFS	☐ 12	$928	☐ 14	$928	☐ 16	$928	☐ 18	$928	☐ 20	$928
		J3PDMTTxxTH	☐ 12	$928	☐ 14	$928	☐ 16	$928	☐ 18	$928	☐ 20	$928
Upper Thoracic (UT)	18.0" (M)	J3PDUTMxxFS			☐ 14	$928	☐ 16	$928	☐ 18	$928	☐ 20	$928
		J3PDUTMxxTH			☐ 14	$928	☐ 16	$928	☐ 18	$928	☐ 20	$928
Shoulder Height (SH)	21.0" (S)	J3PDSHSxxFS			☐ 14	$928	☐ 16	$928	☐ 18	$928	☐ 20	$928
		J3PDSHSxxTH			☐ 14	$928	☐ 16	$928	☐ 18	$928	☐ 20	$928

JAY J3 Back Whitmyer Headrest

AXYS Mount and Pad

☐	J3HEDX8	E0955/E1033	AXYS Headrest with 8" Plush Pad	$449.00
☐	J3HEDX10	E0955/E1033	AXYS Headrest with 10" Plush Pad	$449.00
☐	J3HEDX14	E0955/E1033	AXYS Headrest with 14" Plush Pad	$464.00
☐	J3HEDXAAPN	E0955/E1033	AXYS with Adj-a-Plush Narrow Pad	$566.00
☐	J3HEDXAAPS	E0955/E1033	AXYS with Adj-a-Plush Std Pad	$566.00

LINX Mount and Pad

☐	J3HEDL8	E0955/E1033	LINX Headrest with 8" Plush Pad	$520.00
☐	J3HEDL10	E0955/E1033	LINX Headrest with 10" Plush Pad	$520.00
☐	J3HEDL14	E0955/E1033	LINX Headrest with 14" Plush Pad	$535.00
☐	J3HEDLAAPN	E0955/E1033	LINX with Adj-a-Plush Narrow Pad	$637.00
☐	J3HEDLAAPS	E0955/E1033	LINX with Adj-a-Plush Std Pad	$637.00

JAY J3 Back Laterals *Swing Away Laterals With Pad - Pair*

☐	JLATSA44PR	E0956/E1034	SwingAway Lat with 4"x4" Pad (PR)	$910.00
☐	JLATSA45PR	E0956/E1034	SwingAway Lat with 4"x5" Pad (PR)	$910.00
☐	JLATSA46PR	E0956/E1034	SwingAway Lat with 4"x6" Pad (PR)	$910.00
☐	JLATSA64PR	E0956/E1034	SwingAway Lat with 6"x4" Pad (PR)	$910.00
☐	JLATSA65PR	E0956/E1034	SwingAway Lat with 6"x5" Pad (PR)	$910.00
☐	JLATSA66PR	E0956/E1034	SwingAway Lat with 6"x6" Pad (PR)	$910.00

© 5/2025 MK-100025 REV Q

JAY® Signature Series Cushions

Note: Not all sizes and configurations shown below, for full offering visit www.sunrisemedical.com.

When ordering mark the part number desired and it will be added to the order as a separate line item .

JAY CUSHIONS

JAY Fusion Cushion

✔	Option #	HCPCS	Definition	Price
☐	JFUSION1416	E2622	JAY Fusion 14" x 16"	$515.00
☐	JFUSION1616	E2622	JAY Fusion 16" x 16"	$515.00
☐	JFUSION1618	E2622	JAY Fusion 16" x 18"	$515.00
☐	JFUSION1620	E2622	JAY Fusion 16" x 20"	$515.00
☐	JFUSION1717	E2622	JAY Fusion 17" x 17"	$515.00
☐	JFUSION1816	E2622	JAY Fusion 18" x 16"	$515.00
☐	JFUSION1818	E2622	JAY Fusion 18" x 18"	$515.00
☐	JFUSION1820	E2622	JAY Fusion 18" x 20"	$515.00
☐	JFUSION2018	E2622	JAY Fusion 20" x 18"	$515.00
☐	JFUSION2020	E2622	JAY Fusion 20" x 20"	$515.00
☐	JFUSION2220	E2623	JAY Fusion 22" x 20"	$612.00

JAY Fusion Cushion w/ Single Valve Air Insert

✔	Option #	HCPCS	Definition	Price
☐	JFUSION1416AS	E2622	JAY Fusion 14" x 16"	$578.00
☐	JFUSION1616AS	E2622	JAY Fusion 16" x 16"	$578.00
☐	JFUSION1618AS	E2622	JAY Fusion 16" x 18"	$578.00
☐	JFUSION1620AS	E2622	JAY Fusion 16" x 20"	$578.00
☐	JFUSION1717AS	E2622	JAY Fusion 17" x 17"	$578.00
☐	JFUSION1816AS	E2622	JAY Fusion 18" x 16"	$578.00
☐	JFUSION1818AS	E2622	JAY Fusion 18" x 18"	$578.00
☐	JFUSION1820AS	E2622	JAY Fusion 18" x 20"	$578.00
☐	JFUSION2018AS	E2622	JAY Fusion 20" x 18"	$578.00
☐	JFUSION2020AS	E2622	JAY Fusion 20" x 20"	$578.00
☐	JFUSION2220AS	E2623	JAY Fusion 22" x 20"	$674.00

JAY Fusion Cushion w/ Dual Valve Air Insert

✔	Option #	HCPCS	Definition	Price
☐	JFUSION1416AD	E2622	JAY Fusion 14" x 16"	$583.00
☐	JFUSION1616AD	E2622	JAY Fusion 16" x 16"	$583.00
☐	JFUSION1618AD	E2622	JAY Fusion 16" x 18"	$583.00
☐	JFUSION1620AD	E2622	JAY Fusion 16" x 20"	$583.00
☐	JFUSION1717AD	E2622	JAY Fusion 17" x 17"	$583.00
☐	JFUSION1816AD	E2622	JAY Fusion 18" x 16"	$583.00
☐	JFUSION1818AD	E2622	JAY Fusion 18" x 18"	$583.00
☐	JFUSION1820AD	E2622	JAY Fusion 18" x 20"	$583.00
☐	JFUSION2018AD	E2622	JAY Fusion 20" x 18"	$583.00
☐	JFUSION2020AD	E2622	JAY Fusion 20" x 20"	$583.00
☐	JFUSION2220AD	E2623	JAY Fusion 22" x 20"	$680.00

JAY Union Cushion

✔	Option #	HCPCS	Definition	Price
☐	JU1416	E2607	JAY UNION 14" x 16"	$520.00
☐	JU1616	E2607	JAY UNION 16" x 16"	$520.00
☐	JU1618	E2607	JAY UNION 16" x 18"	$520.00
☐	JU1620	E2607	JAY UNION 16" x 20"	$520.00
☐	JU1717	E2607	JAY UNION 17" x 17"	$520.00
☐	JU1816	E2607	JAY UNION 18" x 16"	$520.00
☐	JU1818	E2607	JAY UNION 18" x 18"	$520.00
☐	JU1820	E2607	JAY UNION 18" x 20"	$520.00
☐	JU2018	E2607	JAY UNION 20" x 18"	$520.00
☐	JU2020	E2607	JAY UNION 20" x 20"	$520.00
☐	JU2220	E2608	JAY UNION 22" x 20"	$635.00

JAY Union Reduced Profile Cushion

✔	Option #	HCPCS	Definition	Price
☐	JU1416RP	E2607	JAY Union RP 14" x 16"	$520.00
☐	JU1616RP	E2607	JAY Union RP 16" x 16"	$520.00
☐	JU1618RP	E2607	JAY Union RP 16" x 18"	$520.00
☐	JU1620RP	E2607	JAY Union RP 16" x 20"	$520.00
☐	JU1717RP	E2607	JAY Union RP 17" x 17"	$520.00

JAY J2 Cushion

✔	Option #	HCPCS	Definition	Price
☐	2108-MJ	E2622	JAY J2 14" x 16"	$524.00
☐	2105-MJ	E2622	JAY J2 15.5" x 16"	$524.00
☐	2106-MJ	E2622	JAY J2 15.5" x 18"	$524.00
☐	2107-MJ	E2622	JAY J2 15.5" x 20"	$553.00
☐	2117-MJ	E2622	JAY J2 17" x 17"	$524.00
☐	2100-MJ	E2622	JAY J2 18" x 16"	$524.00
☐	2101	E2622	JAY J2 18" x 18"	$524.00
☐	2103-MJ	E2622	JAY J2 18" x 20"	$554.00
☐	2118-MJ	E2622	JAY J2 20" x 18"	$554.00
☐	2102-MJ	E2622	JAY J2 20" x 20"	$641.00

JAY J2 Deep Contour Cushion

✔	Option #	HCPCS	Definition	Price
☐	2446	E2622	JAY J2 Deep Contour 14" x 16"	$537.00
☐	2466	E2622	JAY J2 Deep Contour 16" x 16"	$537.00
☐	2468	E2622	JAY J2 Deep Contour 16" x 18"	$537.00
☐	2460	E2622	JAY J2 Deep Contour 16" x 20"	$549.00
☐	2477	E2622	JAY J2 Deep Contour 17" x 17"	$537.00
☐	2486	E2622	JAY J2 Deep Contour 18" x 16"	$537.00
☐	2488	E2622	JAY J2 Deep Contour 18" x 18"	$537.00
☐	2480	E2622	JAY J2 Deep Contour 18" x 20"	$561.00
☐	2408	E2622	JAY J2 Deep Contour 20" x 18"	$561.00
☐	2400-MJ	E2622	JAY J2 Deep Contour 20" x 20"	$677.00
☐	2428	E2623	JAY J2 Deep Contour 22" x 18"	$677.00

JAY X2 Cushion

✔	Option #	HCPCS	Definition	Price
☐	JX21416	E2622	JAY X2 14" x 16"	$501.00
☐	JX21616	E2622	JAY X2 16" x 16"	$501.00
☐	JX21618	E2622	JAY X2 16" x 18"	$501.00
☐	JX21620	E2622	JAY X2 16" x 20"	$501.00
☐	JX21717	E2622	JAY X2 17" x 17"	$501.00
☐	JX21816	E2622	JAY X2 18" x 16"	$501.00
☐	JX21818	E2622	JAY X2 18" x 18"	$501.00
☐	JX21820	E2622	JAY X2 18" x 20"	$501.00
☐	JX22020	E2622	JAY X2 20" x 20"	$501.00

JAY Ion Cushion

✔	Option #	HCPCS	Definition	Price
☐	JION1416	E2607	JAY ION 14" x 16"	$462.00
☐	JION1616	E2607	JAY ION 16" x 16"	$462.00
☐	JION1618	E2607	JAY ION 16" x 18"	$462.00
☐	JION1620	E2607	JAY ION 16" x 20"	$462.00
☐	JION1717	E2607	JAY ION 17" x 17"	$462.00
☐	JION1816	E2607	JAY ION 18" x 16"	$462.00
☐	JION1818	E2607	JAY ION 18" x 18"	$462.00
☐	JION1820	E2607	JAY ION 18" x 20"	$462.00
☐	JION2018	E2607	JAY ION 20" x 18"	$462.00
☐	JION2020	E2607	JAY ION 20" x 20"	$462.00
☐	JION2220	E2608	JAY ION 22" x 20"	$578.00

✔	Option #	HCPCS	Definition	Price
☐	JU1816RP	E2607	JAY Union RP 18" x 16"	$520.00
☐	JU1818RP	E2607	JAY Union RP 18" x 18"	$520.00
☐	JU1820RP	E2607	JAY Union RP 18" x 20"	$520.00
☐	JU2018RP	E2607	JAY Union RP 20" x 18"	$520.00
☐	JU2020RP	E2607	JAY Union RP 20" x 20"	$520.00

Sunrise Medical (US) LLC

12002 Volunteer Blvd. Mount Juliet, TN 37122

Customer Service: 800-333-4000 Fax: 800-300-7502 Email: orders@sunmed.com www.sunrisemedical.com

© 5/2025 MK-100025 REV Q